EN	Evolution-TL and Evolution-C Removal Kits - Instructions for Use	SV	Evolution-TL och Evolution-C borttagningssatser - Bruksanvisning
FR	Trousse d'extraction Evolution-TL et Evolution-C - Mode d'emploi	IS	Evolution-TL og Evolution-C flutningsbúnaður - Notkunarleiðbeiningar
DE	Evolution-TL und Evolution-C Entnahmekits - Gebrauchsanweisung	BG	Комплекти за разглобяване Evolution-TL и Evolution-C Инструкции за употреба
ES	Kits de extracción Evolution-TL y Evolution-C - Instrucciones de uso	HR	Evolution-TL i Evolution-C kompleti za odstranjenje - Upute za korištenje
IT	Kit di rimozione Evolution-TL ed Evolution-C - Istruzioni per l'uso	SR	Evolution-TL и Evolution-C комплекти за одстранивање Упутства за коришћење
PT	Kits de Remoção Evolution-TL e Evolution-C - Instruções de utilização	SL	Evolution-TL in Evolution-C kompleta za odstranjevanje - Navodila za uporabo
CS	Soupravy Evolution-TL a Evolution-C Removal Kits - Návod k použití	RO	Eruse de extracție Evolution-TL și Evolution-C - Instrucțiuni de utilizare
DA	Evolution-TL og Evolution-C Removal Kits - (Fjernelseskits) Brugsanvisning	RU	Наборы экстракторов Evolution-TL и Evolution-C - Инструкция по применению
NL	Evolution-TL- en Evolution-C-verwijderkits - Gebruiksaanwijzing	UK	Набори екстракторів Evolution-TL і Evolution-C - Інструкція по застосуванню
ET	Evolution-TL and Evolution-C Removal Kits - Kasutusjuhend	KK	Evolution -TL мен Evolution-C тасымалды жинақтар Қолданыс туралы нұсқау
FI	Evolution-TL- ja Evolution-C-irrotussarjat - Käyttöohjeet	TK	Evolution-TL ve Evolution-C Çıhaz Çıkartma Kitleri - Kullanım Talimatı
EL	Κιτ αποίεσης Evolution-TL και Evolution-C - Οδηγίες χρήσης	AR	Evolution-C و Evolution أطقم الإزالة تعليمات الاستخدام
HU	Evolution-TL and Evolution-C Removal Kits - Instructions for Use	TC	Evolution-TL Evolution-C 取出套件 使用說明書
LV	Evolution-TL un Evolution-C izņemšanas komplekti - Lietošanas instrukcija	SC	Evolution-TL и Evolution-C 取出套件 使用說明書
LT	„Evolution-TL“ ir „Evolution-C“ stuburo implantų pašalinimo rinkiniai - Naudojimo instrukcijos	KO	Evolution-TL 및 Evolution-C 제거 키트 사용 설명서
NO	Evolution-TL un Evolution-C izņemšanas komplekti - Lietošanas instrukcija	ID	Evolution-TL dan Evolution-C - Petunjuk Penggunaan
PL	Zestawy Evolution-TL Kit oraz Evolution-C Removal Kit-Instrukcja obsługi	JP	Evolution-TL/Evolution-C 抜去キット 取扱説明書
SK	Otočné TL a otočné C extrakčné súpravy Návod na použitie		

Evolution-TL and Evolution-C Removal Kits Instructions for Use

Manufacturer

Tecomet
5307 95th Avenue
Kenosha, WI 53144 USA

Description

The reusable manual orthopedic surgical instruments contained in the Evolution Removal Kits consist of drivers for standard metric and English hex drivers, hexalobe drivers, broken and stripped screw removers, specific drivers, extensions, ratchet handle, multi-axial handle, counter torque wrench, combination wrenches, needle-nose locking pliers and cutting instruments. The instruments are provided in an instrument storage container.

Intended Use

The Evolution-TL Removal Kit is intended for use in the removal of spinal hardware from the lumbar and thoracic regions of the spine. The Evolution-C Removal Kit is intended for use in the removal of spinal hardware from the cervical region of the spine.

The intended use of the various drivers is to remove implanted screws based on the screw head configuration.

The intended use of the ratcheting and fixed drivers is to remove implanted screws, utilizing the appropriate driver types for the implanted screws.

The intended use of the rod cutter is to shear various sized rods based on the clinician and surgical requirements.

The intended use of the open ended wrenches is to remove fasteners.

The intended use of the adjustable pliers is to securely grip and manipulate rods.

The intended use of the extensions is to provide additional length to the various drivers used for the removal of implanted screws.

The intended use of the counter torque wrenches is to provide reactionary torque (counter torque) to a rod during the removal of spinal hardware.

The intended use of the tri-lobe adaptor is to provide a tri-lobe connection to the various drivers used for the removal of implanted screws.

Rx Only

Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare practitioner.

Contraindications

These instruments are contraindicated in the removal of spinal hardware systems which use a pressure fit to retain its components or which twist off the threads of the locking screws.

Precautions

Instruments are provided non-sterile. Clean and sterilize prior to each use.

Neutral pH cleaners are recommended.

Prior to use, care must be taken to protect nerves, vessels and/or organs from damage that may result from the use of these instruments.

Reprocessing Limitations

Repeated processing has minimal effect upon manual instruments. End of life is determined by wear and damage due to use.

Point of Use

Remove excess bodily fluids and tissue with a disposable wipe.

Instruments should not be allowed to dry prior to cleaning.

Containment / Transportation

Follow hospital protocols when handling contaminated and or biohazardous materials. Instrument should be cleaned within 30 minutes after use to minimize the potential of staining, damage, and drying after use.

Manual Cleaning Procedure

1. Prepare a proteolytic enzyme solution according to the manufacturer's instructions.
2. Completely submerge instruments in the enzyme solution and gently shake them to remove trapped bubbles. Actuate instruments with hinges or moving parts to ensure contact of the solution with all surfaces. Lumens, blind holes and cannulations should be flushed with a syringe to remove bubbles and ensure contact of the solution with all instrument surfaces.
3. Soak instruments for a minimum of 10 minutes. While soaking, scrub surfaces using a nylon-bristled brush until all visible soil has been removed. Actuate moveable mechanisms. Particular attention should be given to crevices, hinged joints, box locks, instrument teeth, rough surfaces, and areas with moving components or springs. Lumens, blind holes, and cannulas should be cleaned using a snug fitting round nylon bristle brush. Insert the snug fitting round brush into the lumen, blind hole, or cannula with a twisting motion while pushing in and out multiple times.
NOTE: All scrubbing should be performed below the surface of the enzyme solution to minimize the potential of aerosolizing contaminated solution.
4. Remove the instruments from the enzyme solution and rinse in tap water for a minimum of one (1) minute. Actuate all moveable and hinged parts while rinsing. Thoroughly and aggressively flush lumens, holes, cannulas, and other difficult to access areas.
5. Prepare an ultrasonic cleaning bath with detergent and de-gas according to the manufacturer's recommendations. Completely submerge instruments in the cleaning solution and gently shake them to remove any trapped bubbles. Lumens, blind holes, and cannulations should be flushed with a syringe to remove bubbles and ensure contact of the solution with all instrument surfaces. Sonically clean the instruments at the time, temperature, and frequency recommended by the equipment manufacturer

and optimal for the detergent used. A minimum of ten (10) minutes is recommended.

NOTES:

- **Separate stainless steel instruments from other metal instruments during ultrasonic cleaning to avoid electrolysis.**
 - **Fully open hinged instruments.**
 - **Use wire mesh baskets or tray designed for ultrasonic cleaners.**
 - **Regular monitoring of sonic cleaning performance by means of an ultrasonic activity detector, aluminum foil test, TOSI™ or SonoCheck™ is recommended.**
6. Remove the instruments from the ultrasonic bath and rinse in purified water for a minimum of one (1) minute or until there is no sign of residue detergent or biologic soil. Actuate all moveable and hinged parts while rinsing. Thoroughly and aggressively flush lumens, holes, cannulas, and other difficult to access areas.
 7. Dry instruments with a clean, absorbent, non-shedding lint free cloth. Clean, filtered compressed air may be used to remove moisture from lumens, holes, cannulas, and difficult to access areas.

Automated Cleaning Procedure

A validation study for Automatic Cleaning has not been performed for these instruments. Please consult an appropriate responsible person or hospital protocol for direction.

Disinfection

Disinfection is only acceptable as an adjunct to full sterilization for reusable surgical instruments. See sterilization section below.

Inspection / Functional Testing

1. Carefully inspect each device to ensure that all visible blood and soil has been removed.
2. Visually inspect for damage and/or wear.
3. Check the action of moving parts (such as hinges and box-locks) to ensure smooth operation throughout the intended range of motion.
4. Check instruments with long slender features (particularly rotating instruments) for distortion.
5. Where instruments form part of a larger assembly, check that the devices assemble readily with mating components.
6. Inspect edges of the case, tray and cassette to ensure there aren't any sharp edges which might cut sterilization wrap. Inspect lid to assure there isn't any cracking, bending or warping. Verify lid continues to fit case tightly.

***Note: If damage or wear is noted that may compromise the function of the instrument, do not use the instrument and notify the appropriate responsible person.**

Sterilization

Instruments and the instrument kit must be cleaned prior to sterilization. Instruments must be arranged per the specified tray configuration. Additional instruments not included in the specified tray configuration shall not be added to the tray for sterilization.

With instruments secured within, the kit must be wrapped using CSR (sterilization wrap) following the AAMI double wrap method (ANSI/AAMI ST46), or equivalent. The trays and

lids of the kit must be properly assembled prior to sterilization. It is not recommended to stack assembled kits during sterilization.

Do not place kit on side or upside down. The use of nonabsorbent tray liners can cause condensation to pool.

Method – Pre-Vacuum Steam Sterilization – for 2nd Generation Evolution Systems (Part Numbers ES3XXXX and ES4XXXX ONLY. For example ES30001 or ES40001)

Set parameters at 4 minutes exposure at a temperature of 132°C (270°F) with a minimum of four (4) vacuum pulses prior to actual exposure. A dry time of 30 minutes is recommended.

Method – Gravity Displacement Steam Sterilization – for 2nd Generation Evolution Systems (Part Numbers ES3XXXX or ES4XXXX ONLY. For example ES30001 or ES40001)

Expose instruments for 15 minutes at a temperature of 132° C (270° F). A dry time of 30 minutes is recommended.

Method – Pre-Vacuum Steam Sterilization-for 1st Generation Evolution Systems (Any other Part Number not listed above.)

Set parameters at 10 minutes exposure at a temperature of 132°C (270°F) with a minimum of four (4) vacuum pulses prior to actual exposure.

Method – Gravity Displacement Steam Sterilization- for 1st Generation Evolution Systems (Any other Part Number not listed above)

Expose instruments for 30 minutes at a temperature of 132° C (270° F).

Handle sterilized kit following normal operating room procedures.

Allow case to cool thoroughly for a minimum of 15 minutes after sterilization cycle is complete.

Maintenance

Lubricate hinges, threads and other moving parts with a commercial water-based surgical grade instrument lubricant (such as instrument milk) to reduce friction and wear. Follow lubricant manufacturer's instructions.

Storage

Store sterile, packaged instruments in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, vermin, and temperature/humidity extremes. Examine sterile instrument packages closely prior to opening to ensure that there has been no loss of package integrity.

Instructions for use

Cautions: The user of this instrument set acknowledges that implant manufacturers are constantly improving their designs, which could lead to a fit concern with custom drivers.

Do not use the driver(s) if the instrument does not properly fit into the implanted spinal hardware.

Prior to inserting the driver into the implant, verify that all foreign material is removed from the driving surface. Failure to do so could result in the damaging of implants and/or removal instruments.

The user of the instrument set should ensure that the instrumentation has no more than 60 minutes of direct exposure to the patient. The user of this set should ensure that no

instrumentation/parts of instrumentation be left in the patient during and following the surgical procedure.

Inspect instruments for damage or wear prior to use. Failure to do so could result in the damaging of implants or removal instruments.

***Note: If damage or wear is noted that may compromise the function of the instrument, do not use the instrument and notify the appropriate responsible person.**

The following tables identify the torque ratings of the male hex and male hexalobe drivers. The drivers should not be subjected to torque loads higher than the identified ratings in order to ensure safe use.

Male Hex Driver Torque Ratings


Size	Torque Rating (in.-lbs)
1.5 mm	7
2.0 mm	10
2.5 mm	29
7/64"	29
3.0 mm	47
3.2 mm (1/8")	47
3.5 mm	71
5/32"	100
4.5 mm	138
3/16"	138
5.0 mm	224
7/32"	224
1/4"	224
7.0 mm	224

Male Hexalobe Driver Torque Ratings

Size	Torque Rating (in.-lbs)
T08	10
T10	17
T15	29
T20	47

Size	Torque Rating (in.-lbs)
T25	71
T27	100
T30	138
T40	240


BROKEN SCREW REMOVAL (Combined trephine and Bone Screw Remover)

1. Measure the outside diameter of broken screw and select trephine/ Broken Screw Remover that matches screw size.
2. Place trephine in handle. Turn counter-clockwise around broken screw. The trephine should clear around the screw as well as bite into screw.
3. Continue turning in counter-clockwise direction until screw is completely removed.

Note: *There may be bone in-growth on the screw preventing removal of the screw. Depending on the surgical scenario and condition of implanted screws, the broken screw remover may become damaged after normal usage. Inspect the broken screw remover following usage to determine if any damage has occurred. Discontinue use if damage is present.*

BROKEN SCREW REMOVAL (Separate trephine and Bone Screw Remover)

1. Measure the outside diameter of the screw and select trephine that matches screw size.
2. Place trephine in handle of choice provided in the removal kit. Position trephine over broken screw. Note: The Trephine should clear around the screw. Approximately 5mm of screw should be exposed to properly seat the Trephine. To expose the broken screw, use a clockwise and counter-clockwise motion while applying downward pressure. Advance the trephine in this manner until the trephine reaches the broken screw. Note: During this process the teeth of the trephine may become clogged with debris. Remove trephine and clean debris from tool.
3. Select the bone screw remover that matches the diameter of the bone screw to be removed. Place bone screw remover into handle and position bone screw remover over the exposed screw. Use a counter-clockwise motion while applying downward pressure. The bone screw remover will bite into the screw. Continue turning in a counter-clockwise direction until the screw is completely removed.

Note: There may be bone in-growth on the screw preventing removal of the screw. Depending on the surgical scenario and condition of implanted screws, the broken screw remover may become damaged after normal usage. Inspect the broken screw remover following usage to determine if any damage has occurred. Discontinue use if damage is present.

The male hex drivers, female hex drivers, male hexalobe drivers, and custom drivers are to be used with the appropriate ¼" square or AO ratcheting or fixed driver to remove implanted screws. The driver type should be selected based on the appropriate connection style of the implanted screw. Insert the appropriate driver tip into the mating features of the implanted screw and rotate the driver in the counter-clockwise direction to remove the screw.

The screw extractors can be used to remove stripped implanted screws. The screw extractors are to be used with the appropriate ¼" square or AO ratcheting or fixed driver.

To remove a stripped screw, rotate the screw extractor in the counter-clockwise direction on the stripped screw head surface.

The ratcheting driver is to be used with the appropriate bit driver needed for screw removal. The bit driver can be inserted into the ratcheting driver by compressing the adaptor collar and inserting the appropriate ¼" square or AO driver. Ensure that the driver is fully engaged prior to use. The ratchet can operate in a forward or reverse ratcheting direction, or in a locked setting (fixed driver). The direction of the ratcheting driver can be changed by rotating the selector ring in the clockwise or counter-clockwise direction. The bit driver can be removed from the ratcheting driver by compressing the adaptor collar.

The fixed driver is to be used with the appropriate bit driver needed for screw removal. The bit driver can be inserted into the fixed driver by compressing the adaptor collar and inserting the appropriate ¼" square or AO driver. The bit driver can be removed from the fixed driver by compressing the adaptor collar.

The extensions are to be used with a fixed or ratcheting handle to provide additional length to the bit drivers used to remove implanted screws. The extensions can be inserted into the adaptors of the ratcheting or fixed drivers. Bit drivers can be inserted into the extensions by either compressing the adaptor collar and inserting (AO extensions) or by inserting the male bit driver connection end into the female connection end of the extension (1/4" sq extensions). The bit drivers can be removed from the extensions by either compressing the adaptor collar (AO extensions) or pulling the bit driver out of engagement (1/4" sq extensions).

The counter torque wrenches are used to apply a counter torque or opposite torque to a spinal rod during surgery. The counter torque can be inserted into the appropriate fixed handle. The tip of the counter torque with the full radius slot can be placed over the rod to provide stability as needed.

The adjustable pliers can be used to grip spinal rods during the screw extraction process. The pliers can lock based on the position of the screw knob on one of the handles.

The open ended wrenches are used to remove fasteners. The open ended wrenches can remove various sized fasteners based on the open ended wrench size. Engage the wrench onto the mating features of the implanted screw and rotate the wrench in the counter-clockwise direction to remove the screw.

The tri-lobe adapter can be used with a handle with a female tri-lobe adapter or a 3-jaw chuck adapter. An appropriate sized bit driver can be inserted into the female square end of the tri-lobe adapter in order to remove an implanted screw.

The rod cutter can be used to shear rods in various sizes, as dictated by the surgical procedure and surgeon. To cut a rod, open the rod bender handles by pulling away from each other, insert a rod to be lined up with the rod cutter teeth, and apply force to close the rod cutter handles by pushing the handles towards each other.

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GmbH & Co. KG