IMPORTANT INFORMATION - PLEASE READ BEFORE USE
INFORMATIONS IMPORTANTES - À LIRE AVANT UTILISATION
INFORMAZIONI IMPORTANTI - LEGGERE PRIMA DELL'USO
WICHTIGE INFORMATIONEN - VOR GEBRAUCH BITTE LESEN
INFORMACIÓN IMPORTANTE - LÉALA ANTES DE UTILIZAR ESTOS PRODUCTOS
重要な情報 - 使用前にお読みください
ВАЖНА ИНФОРМАЦИЯ - МОЛЯ, ПРОЧЕТЕТЕ ПРЕДИ УПОТРЕБА
重要信息-使用前请阅读
VAŽNE INFORMACIJE - PROČITATI PRIJE UPORABE
DŮLEŽITÉ INFORMACE - PŘEČTĚTE SI PŘED POUŽITÍM
VIGTIGE OPLYSNINGER - LÆS VENLIGST INDEN BRUG
BELANGRIJKE INFORMATIE - VOOR GEBRUIK DOORLEZEN
OLULINE TEAVE - PALUN LUGEGE ENNE KASUTAMIST
TÄRKEITÄ TIETOJA - LUETTAVA ENNEN KÄYTTÖÄ
ΣΗΜΑΝΤΙΚΕΣ ΠΛΗΡΟΦΟΡΙΕΣ - ΠΑΡΑΚΑΛΟΥΜΕ ΔΙΑΒΑΣΤΕ ΠΡΙΝ ΑΠΟ ΤΗ ΧΡΗΣΗ
FONTOS TUDNIVALÓK - HASZNÁLAT ELŐTT OLVASSA EL
INFORMASI PENTING - HARAP DIBACA SEBELUM MENGGUNAKAN
МАҢЫЗДЫ АҚПАРАТ - ПАЙДАЛАНУ АЛДЫНДА ОҚЫП ШЫҒЫҢЫЗ
중요 정보 - 사용 전에 읽어 보십시오
SVARĪGA INFORMĀCIJA - LŪDZU, IZLASIET PIRMS LIETOŠANAS
SVARBI INFORMACIJA - PRIEŠ NAUDOJANT BŪTINA PERSKAITYTI
VIKTIG INFORMASJON - LES FØR BRUK
WAŻNE INFORMACJE - NALEŻY PRZECZYTAĆ PRZED UŻYCIEM
INFORMAÇÕES IMPORTANTES - LER ANTES DE UTILIZAR
INFORMAȚII IMPORTANTE - VĂ RUGĂM SĂ LE CITIȚI ÎNAINTE DE UTILIZARE
DŮLEŽITÉ INFORMÁCIE - PREČÍTAJTE SI PRED POUŽITÍM
POMEMBNE INFORMACIJE - PREBERITE PRED UPORABO
VIKTIG INFORMATION - LÄS FÖRE ANVÄNDNING
ÖNEMLİ BİLGİ - LÜTFEN KULLANMADAN ÖNCE OKUYUN
ВАЖЛИВА ІНФОРМАЦІЯ - ПРОЧИТАЙТЕ ПЕРЕД ВИКОРИСТАННЯМ!
THÔNG TIN QUAN TRỌNG - VUI LÒNG ĐỌC TRƯỚC KHI SỬ DỤNG

Contents (English - EN).....	1
Contenu (Français - FR)	9
Contenuto (Italiano - IT).....	18
Inhalt (Deutsch - DE).....	27
Índice (Español - ES)	37
内容 (日本語 - JA)	46
Съдържание (български - BG).....	54
内容 (简体中文 - ZH-CN)	64
Sadržaj (Hrvatski - HR)	70
Obsah (Česky - CS).....	78
Indhold (Dansk - DA)	86
Inhoud (Nederlands - NL)	95
Sisukord (Eesti - ET).....	104
Sisältö (Suomi - FI).....	112
Περιεχόμενα (Ελληνικά - EL).....	121
Tartalom (Magyar - HU).....	131
Isi (Bahasa Indonesia - ID).....	140
Мазмұны (Қазақ - KK)	149
내용 (한국어 - KO)	158
Saturs (Latviski - LV).....	165
Turinys (Lietuvių - LT)	173
Innhold (Norsk - NO).....	181
Spis treści (Polski - PL).....	189
Índice (Português - PT)	198
Cuprins (Română - RO)	207
Obsah (Slovenčina - SK).....	216
Vsebina (Slovenščina - SL).....	224
Innehåll (Svenska - SV).....	232
İçindekiler (Türkçe - TR).....	240
Зміст (Українською - UK)	248
Nội dung (Tiếng Việt - VI)	257

Recommended Care, Cleaning, and Sterilization for Reusable Instruments

INSTRUCTIONS FOR USE

These instructions are in accordance with ISO 17664 and apply to:

- Reusable surgical instruments supplied by Tecomet and intended for reprocessing in a health care facility setting. All Tecomet instruments may be safely and effectively reprocessed using the manual or combination manual/automated cleaning instructions provided in this document **UNLESS otherwise noted in instructions accompanying a specific instrument.**

In countries where reprocessing requirements are more stringent than those provided in this document it is the responsibility of the user/processor to comply with those prevailing laws and ordinances.

These reprocessing instructions have been validated as being capable of preparing reusable Tecomet instruments for use. It is the responsibility of the user/hospital/health care provider to ensure that reprocessing is performed using the appropriate equipment and materials, and that personnel have been adequately trained in order to achieve the desired result; this normally requires that equipment and processes are validated and routinely monitored. Any deviation by the user/hospital/health care provider from these instructions should be evaluated for effectiveness to avoid potential adverse consequences.

INSTRUCTIONS FOR USE

Materials & restricted substances

For indication that the device contains a restricted substance or material of animal origin see product label.

Product Description

The Reusable Suction Tips are a family of reusable surgical instruments designed to aid surgeons in rapidly removing bodily and irrigation fluids from a natural or surgical cavity.

Intended Use

These simple manual instruments are attached to a suction system at the distal end of a suction tube and used for the rapid removal of bodily and irrigation fluids from a natural or surgical cavity. They are general-purpose suction tips and are often designed with a finger hole in the handle portion which is used by the operator to control the level of suction applied during the procedure or treatment. They are reusable devices.

Intended Patient Population

The design and intended use of these devices are very generic. They may be used during a wide variety of surgical procedures and are not restricted to a specific surgical patient population.

Indications for Use

The devices are indicated for use when there is a need for rapid removal of bodily and irrigation fluids.

Contra- indications

The devices are for prescription use and should only be used by qualified healthcare personnel. There are no contraindications for the devices.

Intended User

The system or device is prescription-only and therefore should be used by qualified orthopedic surgeons trained in the respective surgical technique.

Accessories and/or other Device(s) intended to be used in combination with the product

For the combinations listed below, ensure firm connection with the assembled device prior to use.

The Reusable Suction Tips are manually connected to sterile suction tubing in the surgical field. The device may be used with a stylet for clearance of obstructions.

Expected Clinical Benefits

When used as intended, the devices assist the surgeon with rapid removal of bodily and irrigation fluids from a natural or surgical cavity.

Adverse Events & Complications

All surgical operations carry risk. The following are frequently encountered adverse events and complications related to having a surgical procedure in general:

- Delay to surgery caused by missing, damaged or worn instruments.
- Tissue injury and additional bone removal due to blunt, damaged or incorrectly positioned instruments.
- Infection and toxicity due to improper processing.

Adverse events to user:

- Cuts, abrasions, contusions or other tissue injury caused by burs, sharp edges, impaction, vibration or jamming of instruments.

Adverse Events & Complications - Reporting of Serious Incidents

Serious Incident Reporting (EU)

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- The death of a patient, user, or other person,
- The temporary or permanent serious deterioration of a patient's, user's, or other person's state of health,
- A serious public health threat.

Where further information is desired, please contact your local Tecomet sales representative. For instruments produced by another legal manufacturer, reference the manufacturer's instructions for use.

Performance and Characteristics

The performance and safety of the system or device have been established, and it represents the current state of the art when used as intended.

Device Life

The life expectancy of reusable suction tips depends on the frequency of use and the care and maintenance they receive. Even with proper handling, care, and maintenance, these devices do not have an indefinite functional life. Therefore, an accurate estimate of their end-of-life is not possible.

Disposal

At the end of the device's life safely dispose of the device in accordance with local procedures and guidelines.

Any device that has been contaminated with potentially infectious substances of human origin (such as bodily fluids) should be handled according to hospital protocol for infectious medical waste. Any disposable or defective device that contains sharp edges should be discarded according to hospital protocol in the appropriate sharps container.

WARNING

- Instruments that are provided NON-STERILE  must be properly cleaned and sterilized prior to each use.
- Ethylene oxide (EO), gas plasma and dry heat sterilization methods are **not recommended** for sterilization of Tecomet reusable instruments. Steam (moist heat) is the recommended method.
- Saline and cleaning/disinfection agents containing aldehyde, chloride, active chlorine, bromine, bromide, iodine or iodide are corrosive and **should not** be used.
- **Do not allow biologic soil to dry on contaminated devices.** All subsequent cleaning and sterilization steps are facilitated by not allowing blood, body fluids and tissue debris to dry on used instruments.
- Automated cleaning using a washer/disinfector alone **may not** be effective for instruments with lumens, blind holes, cannulas, mated surfaces and other complex features. A thorough manual cleaning of such device features is recommended before any automated cleaning process.
- **Metal brushes and scouring pads must not be used** during manual cleaning. These materials will damage the surface and finish of the instruments. Use only soft bristle nylon brushes and pipe cleaners.
- When processing instruments **do not place heavy devices on top of delicate instruments.**
- **Use of hard water should be avoided.** Softened tap water may be used for most rinsing however purified water should be used for final rinsing to prevent mineral deposits.
- Do not process instruments with polymer components at temperatures equal to or greater than 140°C/285°F because severe surface damage to the polymer will occur.
- Oils or silicone lubricants **should not** be used on surgical instruments.
- Prior to cleaning and sterilization, disassemble all instruments that allow for disassembly.
- Personal Protective Equipment (PPE) should be worn when handling or working with contaminated or potentially contaminated instruments.

Precautions

 **ONLY** U.S. Federal law restricts this device to sale by or on the order of a physician.

- The instruments should be inspected for damage and wear before each use. Instruments that show signs of damage or excessive wear should not be used.
- Pay particular attention to quick connections and release buttons, as well as the cutting edges, sharp tips, and teeth.
- Safety caps and other protective packaging material must be removed from the instruments prior to the first cleaning and sterilization.
- As with any surgical instrument, careful attention should be exercised to ensure that excessive force is not placed on the instrument during use. Excessive force can result in instrument failure.
- Careful attention must be paid to asepsis and avoidance of anatomical hazards.
- Caution should be exercised while handling, cleaning, or wiping instruments with sharp cutting edges, tips, and teeth.

Limitations on Reprocessing

- Repeated processing according to these instructions has minimal effect upon metal Tecomet reusable instruments unless otherwise noted. End of life for stainless steel or other metal surgical instruments is generally determined by wear and damage incurred during the intended surgical use.

- Tecomet instruments comprised of polymers or incorporating polymer components can be sterilized using steam however they are not as durable as their metal counterparts. If polymer surfaces show signs of excessive surface damage (e.g. crazing, cracks or delamination), distortion or are visibly warped they should be replaced. Contact your Tecomet representative for your replacement needs.
- Non-foaming, neutral pH enzymatic and cleaning agents are recommended for processing Tecomet reusable instruments. Alkaline agents with a pH of 12 or less may be used to clean stainless steel and polymer instruments in countries where required by law or local ordinance; or where prion diseases such as Transmissible Spongiform Encephalopathy (TSE) and Creutzfeld-Jakob Disease (CJD) are a concern. It is critical that alkaline cleaning agents are completely and thoroughly neutralized and rinsed from the devices or degradation may occur that limits the device life.

REPROCESSING INSTRUCTIONS

Point of Use

- Remove excess biologic soil from the instruments with a disposable wipe. Place devices in a container of distilled water or cover with damp towels.

Note: Soaking in proteolytic enzyme solution prepared according to the manufacturer will facilitate cleaning especially in instruments with complex features such as lumens, mating surfaces, blind holes and cannulas.

- If instruments cannot be soaked or maintained damp then they should be cleaned within 30 minutes of use to minimize the potential for drying prior to cleaning.

Containment and Transportation

- Used instruments must be transported to the central supply for reprocessing in closed or covered containers to prevent unnecessary contamination risk.

Preparation for Cleaning

- Multi-component instruments must be disassembled prior to cleaning. Disassembly, where necessary, is generally self-evident however for more complicated instruments instructions for use are provided and should be followed.
- All cleaning solutions should be prepared at the use-dilution and temperature recommended by the manufacturer. Softened tap water may be used to prepare cleaning solutions.

Note: Fresh cleaning solutions should be prepared when existing solutions become grossly contaminated (turbid).

Manual Cleaning Steps

- **Step 1:** Prepare a proteolytic enzyme solution according to the manufacturer's instructions.
- **Step 2:** Completely submerge instruments in the enzyme solution and gently shake them to remove trapped bubbles. Actuate instruments with hinges or moving parts to ensure contact of the solution with all surfaces. Lumens, blind holes and cannulations should be flushed with a syringe to remove bubbles and ensure contact of the solution with all instrument surfaces.
- **Step 3:** Soak instruments for a minimum of 10 minutes. Scrub surfaces using a soft nylon-bristled brush until all visible soil has been removed. Particular attention should be given to crevices, roughened surfaces, cutting features, hinged joints, sharp edges, box locks, and areas with small components or springs. Lumens, blind holes and cannulas should be cleaned using a long narrow nylon bristle brush/pipe cleaner. Insert a snug fitting long narrow brush/pipe cleaner into the lumen, blind hole or cannula with a twisting motion while pushing in and out multiple times.

- **Step 4:** While immersed, actuate and maneuver all moving parts of the instruments to ensure complete exposure of the cleaning solution to all surfaces.

Note: All scrubbing should be performed below the surface of the enzyme solution to minimize the potential of aerosolizing contaminated solution.

- **Step 5:** Remove the instruments from the enzyme solution and rinse in tap water for a minimum of one (1) minute. Actuate all moveable and hinged parts while rinsing. Thoroughly and aggressively flush lumens, holes, cannulas and other difficult to access areas.
- **Step 6:** Prepare an ultrasonic cleaning bath with detergent and de-gas according to the manufacturer's recommendations. Completely submerge instruments in the cleaning solution and gently shake them to remove any trapped bubbles. Lumens, blind holes and cannulations should be flushed with a syringe to remove bubbles and ensure contact of the solution with all instrument surfaces. Sonically clean the instruments at the time, temperature and frequency recommended by the equipment manufacturer and optimal for the detergent used. A minimum of ten (10) minutes is recommended.

Note: Separate stainless steel instruments from other metal instruments during ultrasonic cleaning to avoid electrolysis. Fully open hinged instruments and use wire mesh baskets or trays designed for ultrasonic cleaners. Regular monitoring of sonic cleaning performance by means of an ultrasonic activity detector, aluminum foil test, TOSI™ or SonoCheck™ is recommended.

- **Step 7:** Remove the instruments from the ultrasonic bath and rinse in purified water for a minimum of one (1) minute or until there is no sign of residue detergent or biologic soil. Actuate all moveable and hinged parts while rinsing. Thoroughly and aggressively flush lumens, holes, cannulas and other difficult to access areas.
- **Step 8:** Visually inspect the instruments for visible soil and repeat the cleaning steps if remaining soil is observed.
- **Step 9:** Dry instruments with a clean, absorbent non-shedding wipe. Clean, filtered compressed air may be used to remove moisture from lumens, holes, cannulas and difficult to access areas.

Combination Manual/Automated Cleaning Steps

- **Step 1:** Prepare a proteolytic enzyme solution according to the manufacturer's instructions.
- **Step 2:** Completely submerge instruments in the enzyme solution and gently shake them to remove trapped bubbles. Actuate instruments with hinges or moving parts to ensure contact of the solution with all surfaces. Lumens, blind holes and cannulations should be flushed with a syringe to remove bubbles and ensure contact of the solution with all instrument surfaces.
- **Step 3:** Soak instruments for a minimum of 10 minutes. Scrub surfaces using a soft nylon-bristled brush until all visible soil has been removed. Particular attention should be given to crevices, hinged joints, roughened surfaces, cutting features, hinged joints, sharp edges, box locks, and areas with small components or springs. Lumens, blind holes and cannulas should be cleaned using a long narrow nylon bristle brush/pipe cleaner. Insert a snug fitting long narrow brush/pipe cleaner into the lumen, blind hole or cannula with a twisting motion while pushing in and out three (3) times.
- **Step 4:** While immersed, actuate and maneuver all moving parts of the instruments to ensure complete exposure of the cleaning solution to all surfaces.

Note: All scrubbing should be performed below the surface of the enzyme solution to minimize the potential of aerosolizing contaminated solution.

- **Step 5:** Remove the instruments from the enzyme solution and rinse in tap water for a minimum of one (1) minute. Actuate all moveable and hinged parts while rinsing. Thoroughly and aggressively flush lumens, holes, cannulas and other difficult to access areas.
- **Step 6:** Place instruments in a suitable validated washer/ disinfectant. Follow the washer/disinfectant manufacturer’s instructions for loading the instruments for maximum cleaning exposure; e.g. open all instruments, place concave instruments on their side or upside down, use baskets and trays designed for washers, place heavier instruments on the bottom of trays and baskets. If the washer/disinfectant is equipped with special racks (e.g. for cannulated instruments) use them according to the manufacturer’s instructions.
- **Step 7:** Process instruments using a standard washer/ disinfectant instrument cycle according to the manufacturer’s instructions. The following minimum wash cycle parameters are recommended:

Cycle	Description
1	Pre-wash ● Cold Softened Tap Water ● 2 minutes
2	Enzyme Spray & Soak ● Hot Softened Tap Water ● 1 minute
3	Rinse ● Cold Softened Tap Water
4	Detergent Wash ● Hot (64-66°C/146-150°F) Tap Water ● 2 minutes
5	Rinse ● Hot (64-66°C/146-150°F) Purified Water ● 1 minute
6	Hot Air Dry (116°C/240°F) ● 7 – 30 minutes

Note: The washer/disinfectant manufacturer’s instructions should be followed. A washer/disinfectant with demonstrated efficacy (e.g. FDA approval, validated to ISO 15883) should be used. Dry time is shown as a range because it is dependent upon the load size placed into the washer/disinfectant. Many manufacturers pre-program their washer/disinfectants with standard cycles and they may include a thermal low-level disinfection rinse after the detergent wash. The thermal rinse may be performed at (82-93°C/180-200°F) and is compatible with Tecomet instruments.

Disinfection

- Tecomet instruments must be terminally sterilized prior to use. See sterilization instructions below.

Drying

- Dry instruments with a clean, absorbent non-shedding wipe. Clean, filtered compressed air may be used to remove moisture from lumens, holes, cannulas and difficult to access areas.

Inspection & Testing

- After cleaning, all devices should be thoroughly inspected for residue biologic soil or detergent. If contamination is still present repeat the cleaning process.
- Visually inspect each device for completeness, damage and excessive wear. If damage or wear is observed that might compromise the function of the device, do not process them further and contact your Tecomet representative for a replacement.
- When inspecting devices look for the following:
 - Cutting edges should be free of nicks and have a continuous edge.
 - Jaws and teeth should align properly.
 - Movable parts should operate smoothly throughout the intended range of motion.
 - Locking mechanisms should fasten securely and close easily.
 - Long thin instruments should be free of bending or distortion.
 - Where instruments form part of a larger assembly, check that all components are available and assemble readily.

Maintenance and Lubrication

- After cleaning and before sterilization, instruments with moving parts (e.g. hinges, box-locks, sliding or rotating parts) should be lubricated with a water-soluble lubricant such as Preserve®, Instrument Milk or equivalent. Always follow the lubricant manufacturer’s instructions for use-dilution, shelf life and application method.

Packaging for Sterilization

- Single devices may be packaged in a medical grade sterilization pouch or wrap. Care should be used when packaging so that the pouch or wrap is not torn. Devices should be wrapped using the AAMI double wrap or equivalent method.
- Reusable wraps are not recommended.
- Instruments may be packaged in rigid general-use trays and cases along with other devices under the following conditions:
 - Arrange all devices to allow access of steam to all surfaces. Open hinged devices and ensure devices are disassembled if they can be.
 - The cases or trays must be wrapped in a medical grade sterilization wrap by following the AAMI double wrap method or equivalent.
 - Total weight of the wrapped case or tray should not exceed 11.4kg/25lbs.
- If the device is sterilized as part of an instrument set in a rigid container, it is the responsibility of the health care facility to ensure that the minimum recommended sterilization parameters are achieved since changes in instrument load size may affect sterilization efficacy. Use only approved sterilization wraps when processing rigid containers that require them.

Sterilization

- Moist heat/steam sterilization is the recommended method for Tecomet instruments.
- Use of an approved chemical indicator within each sterilization load is recommended.
- Always consult and follow the sterilizer manufacturer instructions for load configuration and equipment operation. Sterilizing equipment should have demonstrated efficacy (e.g. FDA approval, process validation, maintenance records).
- Validated exposure times and temperatures to achieve a 10⁻⁶ sterility assurance level (SAL) are listed in the following table.

Cycle Type	Temperature	Exposure Time (Wrapped)
Pre-vacuum / Vacuum Pulse	132°C / 270°F	4 minutes
Pre-vacuum / Vacuum Pulse	134°C / 273°F	3 minutes

Note: A 30 minute minimum cooling time is recommended after drying but longer times may be necessary because of load configuration, ambient temperature and humidity, device design and packaging used.

Note: Drying times may need to be increased for larger loads.

Note: Disinfection/steam sterilization parameters recommended by the World Health Organization (WHO) for reprocessing instruments where there is a concern about TSE/CJD contamination are: 134°C/273°F for 18 minutes.

Storage

- Sterile packaged instruments should be stored in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, vermin and temperature/humidity extremes.

Note: Inspect every package before use to ensure that the sterile barrier (e.g. wrap, pouch, filter) is not torn, perforated, shows signs of moisture or appears to be tampered with. If any of those conditions are present then the contents are considered non-sterile and should be re-processed through cleaning, packaging and sterilization.

SYMBOLS USED ON LABELING

Caution

Non-Sterile

U.S Federal law restricts this device to sale by or on the order of a physician

CE Mark with Notified Body #

Authorized Representative in the European Community

Manufacturer

Date of Manufacture

Lot Number

Catalog Number

Consult Instructions for Use

Medical Device

Packaging Unit

Country of Manufacture

Distributor

Swiss Authorized Representative

Importer

Unique Device Identifier

Additional Labeling

“MATERIAL”

Raw materials utilized in the manufacturing of the device

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Symmetry Medical Manufacturing, Inc.
Medtronic Sofamor Danek USA, Inc.
Symmetry Medical Polyvac S.A.S