IMPORTANT INFORMATION - PLEASE READ BEFORE USE  
INFORMATIONS IMPORTANTES - À LIRE AVANT UTILISATION  
INFORMAZIONI IMPORTANTI - LEGGERE PRIMA DELL'USO  
WICHTIGE INFORMATIONEN - VOR GEBRAUCH BITTE LESEN  
INFORMACIÓN IMPORTANTE - LÉALA ANTES DE UTILIZAR ESTE PRODUCTO  
重要な情報 - 使用前にお読みください  
ВАЖНА ИНФОРМАЦИЯ - МОЛЯ, ПРОЧЕТЕТЕ ПРЕДИ УПОТРЕБА  
VAŽNE INFORMACIJE - PROČITATI PRIJE UPORABE  
DŮLEŽITÉ INFORMACE - PŘEČTĚTE SI PŘED POUŽITÍM  
VIGTIGE OPLYSNINGER - SKAL LÆSES INDEN BRUG  
BELANGRIJKE INFORMATIE - VOOR GEBRUIK DOORLEZEN  
OLULINE TEAVE - PALUN LUGEGE ENNE KASUTAMIST  
TÄRKEITÄ TIETOJA - LUETTAVA ENNEN KÄYTTÖÄ  
ΣΗΜΑΝΤΙΚΕΣ ΠΛΗΡΟΦΟΡΙΕΣ - ΠΑΡΑΚΑΛΟΥΜΕ ΔΙΑΒΑΣΤΕ ΠΡΙΝ ΑΠΟ ΤΗ ΧΡΗΣΗ  
FONTOS TUDNIVALÓK - HASZNÁLAT ELŐTT OLVASSA EL  
SVARĪGA INFORMĀCIJA - LŪDZU, IZLASIET PIRMS LIETOŠANAS  
SVARBI INFORMACIJA - PERSKAITYKITE PRIEŠ NAUDODAMI  
VIKTIG INFORMASJON - LES FØR BRUK  
WAŽNE INFORMACJE - NALEŻY PRZECZYTAĆ PRZED UŻYCIEM  
INFORMAÇÃO IMPORTANTE - LER ANTES DE UTILIZAR  
INFORMAȚII IMPORTANTE - VĂ RUGĂM SĂ LE CITIȚI ÎNAINTEA UTILIZĂRII  
ВАЖНАЯ ИНФОРМАЦИЯ - ИЗУЧИТЕ ПЕРЕД ПРИМЕНЕНИЕМ ИЗДЕЛИЯ  
VAŽNE INFORMACIJE - PROČITATI PRIJE UPORABE  
DŮLEŽITÉ INFORMÁCIE - PREČÍTAJTE SI ICH PRED POUŽITÍM  
POMEMBNE INFORMACIJE - PREBERITE PRED UPORABO  
VIKTIG INFORMATION - LÄS FÖRE ANVÄNDNING  
ÖNEMLI BİLGİ - LÜTFEN KULLANMADAN ÖNCE OKUYUN

Contents (English - EN).....	1
Contenu (Français - FR).....	5
Contenuto (Italiano - IT).....	9
Inhalt (Deutsch - DE).....	13
Índice (Español - ES).....	17
内容 (日本語 - JA).....	21
Съдържание (Български - BG).....	25
Sadržaj (Hrvatski - HR).....	29
Obsah (Česky - CS).....	33
Indhold (Dansk - DA).....	37
Inhoud (Nederlands - NL).....	41
Sisukord (Eesti - ET).....	45
Sisältö (Suomi - FI).....	49
Περιεχόμενα (Ελληνικά - EL).....	53
Tartalom (Magyar - HU).....	57
Saturs (Latviski - LV).....	61
Turinys (Lietuvių - LT).....	65
Innhold (Norsk - NO).....	69
Spis treści (Polski - PL).....	73
Índice (Português - PT).....	77
Cuprins (Română - RO).....	81
Содержание (Русский - RU).....	85
Sadržaj (Srpski - SR).....	89
Obsah (Slovenčina - SK).....	93
Vsebina (Slovenščina - SL).....	97
Innehåll (Svenska - SV).....	101
İçindekiler (Türkçe - TR).....	105

## **INSTRUCTIONS FOR USE/INTENDED USE**

The X Series Power System Instrument Cases are intended for the protection, organization, and delivery to the surgical field of surgical instruments and/or other medical devices. The X Series Power System Cases are not intended to maintain sterility by itself. It is designed to facilitate the sterilization process when used with an FDA approved sterilization wrap or rigid sterilization container. Wrapping material is designed to allow air removal, steam penetration/evacuation (drying) and to maintain the sterility of the internal components.

### **Recommendations for Care, Cleaning, and Sterilization for the X Series Power System Instruments Case**

Both physical and chemical (detergent) processes may be necessary to clean soiled items. Chemical (detergent) cleaners alone cannot remove all soil and debris; therefore, a careful manual cleaning of each item with a soft sponge or cloth is essential for maximum decontamination. For difficult areas, a clean soft bristled brush is recommended. Once the items have been cleaned, they should be thoroughly rinsed with clean water to remove any detergent or chemical residue before sterilization. **Tecomet** recommends the use of a mild enzymatic detergent with a near neutral pH. Do not use solvents, abrasive cleaners metal brushes, or abrasive pads. Baskets may be placed in mechanical cleaning equipment.

**Always** inspect instruments and case for cleanliness and damage before use. Make sure all latches and handles are secure and in working order. Do not overload baskets. Tecomet recommends that the baskets be wrapped according to the sterilization wrap manufacturer's instructions prior to sterilization to maintain sterility of internal components/items and for proper aseptic presentation to the surgical field. Always follow the X Series Owner's Manual if their sterilization or drying recommendations exceed these guidelines. For steam autoclaves, Tecomet has verified product performance in the following cycles.

## Pre-vacuum Sterilizer

Wrapped baskets should be sterilized according to the table below.

Cycle Type	Temperature	Exposure Time	Dry Time	Cool-Down Time*
<b>United States Recommended Parameters</b>				
Pre-vacuum / Vacuum Pulse	132°C / 270°F	4 minutes	45 minutes	30 minutes
Cycle Type	Temperature	Exposure Time	Dry Time	Cool-Down Time*
<b>European Recommended Parameters</b>				
Pre-vacuum / Vacuum Pulse	134°C / 273°F	3 minutes	45 minutes	30 minutes

\*Allow the wrapped device to cool outside of the autoclave for 30 minutes prior to use.

## Filtered Sterilization Container Systems

Tecomet delivery systems have not been validated for use in filtered sterilization containers. Tecomet does not recommend use of such a system.

If filtered sterilization container systems are used, the user is responsible for following the recommendations of the manufacturer for the proper placement and use of cases and tray inside the container.

## Contraindications

The X Series Power System Instrument Cases have **NOT** been validated for use with flexible endoscopes, or devices with lumens or working channels longer than 4 inches (3mm ID). Always refer to the X Series Owner's Manual. The X Series Power System Instrument Cases have **NOT** been validated for ETO sterilization of devices. Always refer to the X Series Owner's Manual. Stacking the X Series Power System Instrument Cases and overloading the cases will adversely affect sterilization and drying effectiveness. **DO NOT STACK** cases or trays in an autoclave chamber.

## Notice

1. **DO NOT** load cases into sterilizer on sides or upside down with lid side on shelf or cart. Load cases on cart or shelf so that the lid is always facing upwards. This will allow for proper drying. The X Series Power System Instrument Cases are designed to drain in this position.
2. **DO NOT STACK** cases or trays in an autoclave chamber.
3. After the autoclave door is opened, all cases must be allowed to cool thoroughly. Place cases on rack or shelf with linen cover until cooling is complete. The potential for condensation may increase if the case is not allowed to cool properly.
4. If condensation is observed, check to ensure steps 1, 2, and 3 were followed. In addition, verify that the steam, which is being used for sterilization processing, has a quality of more than 97%<sup>1</sup>. Also confirm that the sterilizers have been inspected for routine maintenance in accordance with manufacturer recommendations.

### Symbols Used on Labeling:

Caution

Non-Sterile

U.S. Federal law restricts this device to sale by or on the order of a physician

CE Mark

Authorized Representative in the European Community

Manufacturer

Date of Manufacture

Lot Number

Catalog Number

Consult Instructions for Use

European Union Importer

Medical Device

Unique Device Identifier

Packaging Unit

Country of Manufacture

Distributor

**The following table provides guidance on maximum loads. Do not overload delivery systems (baskets), and always follow AAMI, OSHA and hospital standards for maximum loading. Follow the X Series Owner's Manual's instructions for use.**

Tray Size	Maximum Total Weight
Small (~ 4 x 8 x 1 in.)	1 lb. (0.45 kg)
Tray (~ 8 x 17 x 1 in.)	2 lb. (0.91 kg)
½ Size Case (~ 7 x 7 x 4 in.)	8 lb. (3.6 kg)
¾ Size Case (~ 7 x 12 x 4 in.)	13 lb. (5.9 kg)
Full Size Case (~ 9 x 19 x 4 in.)	22 lb. (10.0 kg)

## REFERENCES

<sup>1</sup>ANSI/AAMI ST79: Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities

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Symmetry Medical Manufacturing, Inc.
Symmetry Medical Polyvac S.A.S