

# CROSSBAR ACETABULAR REAMERS

Our **Cross Bar Acetabular Reamers** are the newest solution for preparing the acetabulum during Total Hip Arthroplasty (THA) surgery.

The cross bar acetabular reamer is the latest addition to Tecomet's all encompassing hip preparation instrument product line and continues Tecomet's innovative reputation with excellent cutting performance and overall longevity.

- Designed to be attached to acetabular reamer drivers accepting a cross bar attachment—this includes both Tecomet designs and competitive drivers
- Case Delivery System available—accommodates 40 reamers
- Reamer tooth configuration is designed to produce a smooth hemispherical acetabular surface
- Tecomet's Acetabular Reamers are not specific to any Acetabular Shell/Cup Implant
- Cleaning & Sterilization methods are outlined in Tecomet's IFU #7033-99
- Structural integrity of the cross bar weld connection to the acetabular reamer dome have been verified per successful static torque test TS-19-0111
- Utilizes an exclusive 420 SS material composition based on previous acetabular reamer design offerings

Contact your local Sales Manager or visit [tecomet.com](http://tecomet.com) for more information.



3724 North SR 15  
Warsaw, IN 46582 USA  
Phone +1 574.267.8700

[sales@tecomet.com](mailto:sales@tecomet.com)  
[www.tecomet.com](http://www.tecomet.com)

Cat #: 5103-00  
Rev. E  
Printed in the USA ©2023  
Tecomet. All rights reserved.

# CROSSBAR ACETABULAR REAMERS

## Part Numbers

36 mm	16-077-036	51 mm	16-077-051	66 mm	16-077-066
37 mm	16-077-037	52 mm	16-077-052	67 mm	16-077-067
38 mm	16-077-038	53 mm	16-077-053	68 mm	16-077-068
39 mm	16-077-039	54 mm	16-077-054	69 mm	16-077-069
40 mm	16-077-040	55 mm	16-077-055	70 mm	16-077-070
41 mm	16-077-041	56 mm	16-077-056	71 mm	16-077-071
42 mm	16-077-042	57 mm	16-077-057	72 mm	16-077-072
43 mm	16-077-043	58 mm	16-077-058	73 mm	16-077-073
44 mm	16-077-044	59 mm	16-077-059	74 mm	16-077-074
45 mm	16-077-045	60 mm	16-077-060	75 mm	16-077-075
46 mm	16-077-046	61 mm	16-077-061	76 mm	16-077-076
47 mm	16-077-047	62 mm	16-077-062	77 mm	16-077-077
48 mm	16-077-048	63 mm	16-077-063	78 mm	16-077-078
49 mm	16-077-049	64 mm	16-077-064	79 mm	16-077-079
50 mm	16-077-050	65 mm	16-077-065	80 mm	16-077-080



3724 North SR 15  
Warsaw, IN 46582 USA  
Phone +1 574.267.8700

[sales@tecomet.com](mailto:sales@tecomet.com)  
[www.tecomet.com](http://www.tecomet.com)

Cat #: 5103-00  
Rev. E  
Printed in the USA ©2023  
Tecomet. All rights reserved.