



EVOLUTION™ – C KIT

Cervical Implant Removal Systems

Second Generation with additional bit/drivers for further hardware removal solutions:

Custom Drivers ♦ Male & Female Hex Drivers ♦ Hexalobes ♦ Broken Screw Extractors



Contact your local Sales Manager or visit tecomet.com for more information.

Evolution™ -C Kit Part #: ES40001

Evolution™ kits are designed to remove/revise a wide range of the major thoracolumbar or cervical implant systems on the market. Our comprehensive implant removal systems contain virtually all the bits and instruments that will optimize the surgeon's removal options—based on both procedural and surgical preference. The compact kit design also features StackPak® trays and lids for easy assembly and transportation.

Top tray efficiently organizes standard metric and custom drivers

Top Tray Part #:
TY20001

Bottom tray organizes the necessary instruments, handles, extensions, and wrenches

Bottom Tray Part #:
TY20002

Improved vent to volume ratio for better steam sterilization

Case, Lid Part #:
LC50001



EVOLUTION™ – C KIT

Part Numbers



Hex Drivers Size (Male)



1.5mm	HD2060591
2.0mm	HD2060787 (2)
2.5mm	HD2060984 (2)
3.0mm	HD2061181
3.2mm	HD2061250
3.5mm	HD2061378



Set Screw Size Extractors



1.5mm	RS1020591
2.0mm	RS1020787
2.5mm	RS1020984
3.0mm	RS1021181
3.5mm	RS1021378
4.0mm	RS1021575



Cruciforms



Ø4.5mm	CD102002
Ø3.0mm	CD102003
Ø6.0mm	CD102011



Trehines Size



3.5mm	HR1041378
4.0mm	HR1041575
4.5mm	HR1041771
5.0mm	HR1041969



Hexalobe Drivers Size



T8	HD207T08
T10	HD207T10
T15	HD207T15
T20	HD207T20
T25	HD207T25









Blade









Flat Blade, 1mm x 4mm	CD102001
--------------------------	----------



Custom Drivers







	Five-Star Driver	CD102004
	Tri-Lobe Driver, A-O	CD102005
	Two-Prong, A-O	CD102006
	Hexalobe, A-O	CD102007-X15 (2)
	Hexalobe, A-O	CD102007-X20
	Hexalobe, A-O	CD102007-X25
	Cam, A-O	CD102008

Custom Drivers, cont.

	Outer Nut, Dodecagon	CD102009
	Quad Driver	CD102012
	Cervive® Three-Prong	CD102013
	Driver for Cervi-Lok® System	CD102014
	Driver for Spinelink® ACS System	CD102015
	Cam Tightener DePuy®	CD102016

Unlocking Tools

	Unlocking tool for Seaspine® System	UT101001
	Unlocking tool for Alphatec® System	CD102010

Instruments		
	Adjustable, Needle Nose Pliers	AP211001
	Shaft Extension	EX201076 (2)
	Axial Mini Handle, A-O	MA101000
	SD, Ratchet, A-O	SR632000
	T-Handle, A-O	TA101000
	4mm Rod Cutter	46629