

HIGH PERFORMANCE ROD CUTTER

The customizable **High Performance Rod Cutter** gives you clean, consistent cuts on cobalt chromium, 316 stainless steel, or titanium with less effort.

Ergonomic handles: Ortho-Grip® silicone handles provide a comfortable, slip-proof grip.

New engineering, better results: Designed with an optimized mechanical advantage to increase leverage and momentum, our proprietary shearing technology and breakthrough head design result in clean cuts every time.

Adjustable arms: Handles extend up to 28" and lock in place for increased leverage and reduced operating effort.

Rod guide for easier alignment: With holes for four different size rods, the guide allows for easy insertion and proper positioning in the cutting hole while minimizing slack.

Contact your local Sales Manager or visit tecomet.com for more information.



5307 95th Avenue
Kenosha, WI 53144
Phone +1 262.925.4500

sales@tecomet.com
www.tecomet.com

Cat #: 5103-00
Rev. E
Printed in the USA ©2023
Tecomet. All rights reserved.

HIGH PERFORMANCE ROD CUTTER

Other Variations



Handheld Rod Cutter

- Dual-lever cutting mechanisms for easy, burr-free cutting of cobalt chrome, titanium or stainless steel rods from 2.50mm to 6.35mm in diameter
- Locking handle extends to provide increased leverage for cutting, and collapses for storage and autoclaving
- Two options: Base mounted or removable base
- Brushed, matte, or glass bead or stainless steel—or any suitable material and finish per customer specifications
- Ergonomic Ortho-Grip® silicone rubber handles for comfortable grip and slip reduction are available in virtually any Pantone® color



Open Ended Cutter

- Easy, burr-free cutting of stainless steel rods via perpendicular approach to the rod. Ideal for cutting stepped rod or custom implants
- Distal tip has an open-ended mechanism design intended to allow for easier access while cutting
- Large profile cutting jaw ensures maximum strength for cutting
- Locking handle extends to provide increased leverage for cutting, and collapses for storage and autoclaving
- Ergonomic Ortho-Grip® silicone rubber handles for comfortable grip and slip reduction are available in virtually any Pantone® color